


# Production and actions of small RNAs

V. Narry Kim and Mikiko C. Siomi

Small (20–30 nt) RNAs are associated with members of the Argonaute (Ago) family, which comprises two subfamilies: Ago and Piwi. Based on their biogenesis mechanism and the type of Argonaute proteins that they associate with, at least three classes of small RNAs can be distinguished in eukaryotes: microRNAs (miRNAs), endogenous small interfering

RNAs (endo-siRNAs) and Piwi-interacting RNAs (piRNAs). miRNAs control mRNA stability and translation by targeting cognate mRNAs. endosiRNAs suppress repetitive genes by cleaving their transcripts. Some piRNAs mediate RNA cleavage or heterochromatin formation of transposons, although the functions of most piRNAs are still unknown.








## Possible mechanisms of action

with Ago-subfamily proteins. In humans, endo-siRNAs are thought to be

loaded onto AGO2; it has not been determined whether they associate with








abundant groups of miRNA-like RNAs that are derived from various sources

of short hairpins in a Drosha-independent manner.



				CCR4-NOT	HP1 HP1
Table   Eukaryotic small RNAs are associated with Argonaute-family proteins					
Subfamily	Ago-family protein	Class of small RNA*	Length of small RNA	Origin of small RNA <sup>‡</sup>	Mechanism of action
Mammals					
Ago	AGO1–4	miRNA	21–23 nt	miRNA genes	Translational repression, mRNA degradation, mRNA cleavage and heterochromatin formation?
		endo-siRNA§	21–22 nt	Intergenic repetitive genes, pseudogenes and endo-siRNA clusters	mRNA cleavage?
Piwi	MILI (PIWIL2 in humans)	Pre-pachytene piRNA and pachytene piRNA	24–28 nt	Transposons and piRNA clusters	Heterochromatin formation (DNA methylation)
	MIWI (PIWIL1 in humans)	Pachytene piRNA	29–31 nt	piRNA clusters	?
	MIWI2 (PIWIL4 in humans)	Pre-pachytene piRNA	27–29 nt	Transposons and piRNA clusters	Heterochromatin formation (DNA methylation)
	(PIWIL3 in humans)	?	?	?	?
Drosophila melanogaster					
Ago	AGO1	miRNA	21–23 nt	miRNA genes	Translational repression and mRNA degradation
	AGO2	_endo-siRNA	~21 nt	Transposons, mRNAs and repeats	RNA cleavage
		exo-siRNA	~21 nt	Viral genome	Viral RNA cleavage
Piwi	AUB	piRNA	23–27 nt	Transposons, repeats, piRNA clusters and Su(Ste) locus	RNA cleavage
	AGO3	piRNA	24–27 nt	Transposons and repeats (unknown in testis)	RNA cleavage
	PIWI	piRNA	24–29 nt	Transposons, repeats and piRNA clusters	Heterochromatin formation?
Schizosaccharomyces pombe Schizosaccharomyces pombe					
Ago	Ago1	endo-siRNA	~21 nt	Outer centromeric repeats, mating-type locus and subtelomeric regions	Heterochromatin formation
Arabidopsis thaliana <sup>  </sup>					
Ago	AGO1	miRNA	20–24 nt	miRNA genes	mRNA cleavage and translational repression
			21 nt	TAS genes	mRNA cleavage
		exo-siRNA	20–22 nt	Viral genome	Viral RNA cleavage
	AGO4 and AGO6	rasiRNA	24 nt	Transposons and repetitive elements	Heterochromatin formation
	AGO7	miR-390	21 nt	miRNA gene	Cleavage of TAS3 RNA

\*Small RNAs that are the main partners of a given Ago protein are listed. †miRNAs, as a class, are expressed in all cell types, whereas endo-siRNAs and piRNAs are expressed abundantly in germline cells and contribute to germline development. §So far, only AGO2 has been shown to be required for endo-siRNAs. Plants have ten Ago proteins, but only those with known small RNA partners are shown.

## Components of the RISC complex at Abcam

Abcam is able to deliver a comprehensive portfolio of all the very best and most up to date antibodies to components of the RISC complex. We are constantly increasing the number of ChIP-grade and batch tested antibodies, enabling the research community to deconstruct epigenetic pathways.

Three of our very best Dicer and Ago antibodies:

• Ago1 (Drosophila) (ab5070)

• DDX6 (ab40684)

- Ago2 (Drosophila) (ab5072)
- Two of our very best RISC antibodies: • TRBP (ab42018)

The Abpromise • Dicer 1 (ab13502)

If our antibodies do not perform as described on the datasheet, notify us within 4 months of delivery so we can help you or offer a replacement or a refund.

All our antibodies to the components of the RISC complex are

highly characterized and our datasheets provide a library of

Investigate these and other chromatin and RNAi antibodies

relevant information.

further at www.abcam.com/RISC



## **Abbreviations**

Ago, Argonaute; AUB, Aubergine; CCR4, C-C chemokine type 4; DGCR8, DiGeorge syndrome critical region gene 8; DNMT, DNA methyltransferase; dsRNA, double-stranded RNA; endo-siRNA; endogenous small interfering RNA; exo-siRNA, exogenous small interfering RNA; H, histone; HMT, histone methyltransferase; HP1, heterochromatin protein 1; LOQS, Loquacious; m<sup>7</sup>G, 7-methylquanosine; Me, methyl; miRNA, microRNA; nt, nucleotide; piRNA, Piwi-interacting RNA; Pol II, RNA polymerase II; PACT, PKR-activating protein; pre-miRNA, precursor miRNA; pri-miRNA, primary miRNA; rasiRNA, repeat-associated small interfering RNA; ssRNA, single-stranded RNA; Su(Ste), Suppressor of Stellate; TAS, tasi gene; tasiRNA, trans-acting siRNA; TRBP, HIV-1 TAR RNA-binding protein.

Poster design by Vicky Askew, edited by Arianne Heinrichs, copyedited by Simon Bishop. © 2009 Nature Publishing Group.

## **Contact information**

V. Narry Kim is at the School of Biological Sciences and Center for National Creative Research, Seoul National University, Seoul, 151-742, South Korea. | Mikiko C. Siomi is at Keio University School of Medicine, 35 Shinanomachi, Shinjuku-ku, Tokyo 160-8582, and at the Japan Science and Technology Agency (JST), Core Research for Evolutional Science and Technology (CREST), Saitama 332-0012, Japan. e-mails: narrykim@snu.ac.kr; siomim@sc.itc.keio.ac.jp

For further reading, see www.nature.com/nrm/posters/smallrnas

## Linked review article

Kim. V. N., Han, I. & Siomi, M. C. Biogenesis of small RNAs in animals. Nature Rev. Mol. Cell Biol. (doi:10.1038/nrm2608).